

Installing Cleverscope

Thanks very much for purchasing a Cleverscope!

Please follow these instructions to install Cleverscope:

1. Run setup.exe on the CD that was supplied with the Cleverscope. Follow the prompts, and Cleverscope should be installed. Leave the CD in the CD drive.
2. Connect up the acquisition unit - you'll need to plug in the power adaptor, and the USB cable to a spare port.
3. The PC should discover the USB connection and display 'Cleverscope CS328 scope' or similar. Just wait after this for Windows to load the driver request dialog. (It might take some seconds as some installations go looking on the internet for a driver). Once the dialog comes up, choose to find the driver automatically. After some searching Windows will find the driver on the CD and install it. It will warn you that it is not digitally signed. Just carry on - we are still waiting for the digital certificate, and no harm will come to your machine. After some more work Windows should announce that the hardware is installed.
4. Startup Cleverscope (it will be in the Start Menu. You can drag it to the desktop if desired).
5. Click on 'Auto' in the Cleverscope Control Panel, and the graph should become live.
6. Exit Cleverscope (File/Exit) or close button.
7. Now follow the section 'Using the CS328 acquisition unit' , page 18 of the manual to set up the probes.

Any problems, please email bart.schroder@cleverscope.com

You may wish to register any suggestions on the forum at <http://www.cleverscope.com/bb/>.

You will find some useful files on the CD. These are:

- Cleverscope Guided Tour.pdf This gives you a quick overview of what you can do
- Cleverscope Screen Shots.pdf Shows you screen shots for common functions
- Cleverscope Manual v2-2.pdf The Cleverscope user manual
- Cleverscope Quick Ref.pdf This is a quick reference sheet for common activities